

WORK INSTRUCTION

TITLE: Product Label Approval Process
--

1.0 PURPOSE AND SCOPE

- 1.1 To allow Maxx Orthopedics the ability to approve labels prior to use on product. All parts being packaged for Maxx Orthopedics will need to have label copy approved by Maxx Orthopedics quality representative **prior** to being placed on finished product.

2.0 REFERENCES

- 2.1 MXO-PP00001 Labeling Guidelines, Freedom Total Knee System & Libertas Hip System
- 2.2 MXO-PP00006 – Maxx Orthopedics GTIN Codes
- 2.3 WI051
- 2.4 WI052-F1

3.0 PROCEDURE FOR LABELS FROM MERIL

- 3.1 Production and Quality Control at Meril Healthcare will perform check of labels to be used on Maxx Orthopedics product prior to use.
- 3.2 Meril creates and prints the labels, printing an extra label for each batch, and the additional label is placed on Meril’s “Reference Label 1st Copy Before Printing Lot Label” form. The label is reviewed by a representative from Meril’s Production team and signed in the “Checked By Production” location.
- 3.3 A representative of Meril’s Quality Control team will review the label and will sign in the “Verified By QC” location prior to sending copy to Maxx Orthopedics.
- 3.4 The “Reference Label 1st Copy Before Printing Lot Label” form will be emailed to Maxx Orthopedics quality team at: maxxquality@maxxortho.net **prior** to label application to finished product. This will help in finding any errors that labels may have and is required as part of Maxx Orthopedics label requirements as listed in section 2.0.

WORK INSTRUCTION

TITLE: Product Label Approval Process
--

- 3.5 A representative of Maxx Orthopedics quality team will review label copy and compare them to corresponding Production Procedures referenced in section 2.0 to ensure all necessary and correct information is on the labels.
- 3.6 Maxx Orthopedics quality representative will also scan the QR codes to ensure the accuracy of information is contained in the QR codes as well.
- 3.7 If approved, the Maxx Orthopedics quality representative will sign in the “Maxx Approval” location of the form.
- 3.8 The Maxx Orthopedics quality representative will then scan the approved “Reference Label 1st Copy Before Printing Lot Label” form and will send the confirmation to Meril via email that the labels are approved and can be affixed to the product/packaging. The “Reference Label 1st Copy Before Printing Lot Label” form will then become part of the DHR (Device History Record).
- 3.9 If not approved, Maxx Orthopedics will communicate to Meril the necessary changes needed to the label. Once changed, Meril will send the corrected labels to Maxx Orthopedics again for approval.

4.0 PROCEDURE FOR LABELS FROM ACI

- 4.1 ACI creates and prints the labels, printing an extra label for each batch, and the additional label is placed on a “Label Approval Record” form (ACI Form # 0014)
- 4.2 Representatives from Manufacturing and Quality Assurance at ACI review the label and sign in the “Approvals” section prior to sending copy to Maxx Orthopedics.
- 4.3 The “Label Approval Record” form will be emailed to the Maxx Orthopedics quality team maxxquality@maxxortho.net **prior** to label application to the finished product. This will help in finding any errors labels may have and is required as part of Maxx Orthopedics label requirements as listed in section 2.0.

WORK INSTRUCTION

TITLE: Product Label Approval Process
--

- 4.4 A representative of the Maxx Orthopedics quality team will review the label copy and compare them to corresponding Production Procedures referenced in section 2.0 to ensure all necessary and correct information is on the labels.
- 4.5 A Maxx Orthopedics quality representative will also scan the QR codes to ensure the accuracy of information is contained in the QR codes as well.
- 4.6 If approved, a Maxx Orthopedics quality representative will email ACI approval to use/affix the labels on the product/packaging.
- 4.7 If not approved, Maxx Orthopedics will communicate to ACI the necessary changes needed to the label. Once changed, ACI will send the corrected labels to Maxx Orthopedics again for approval.

5.0 PROCEDURE FOR LABELS PRINTED AT MAXX ORTHOPEDICS

- 5.1 A representative from the Maxx Orthopedics Production department will create the label and print labels per Work Instruction WI051. An additional label is printed for each label print batch.
- 5.2 The Production representative will place the additional label printed on Maxx Orthopedics Label Approval Record (WI052-F1) and complete the product information at the top of the form:
 - Item # - Part number (REF) on the product label
 - Qty – Quantity of parts to be labeled in the batch
 - Lot # - Lot number (LOT) on the product label
 - Description – Product description per product label
- 5.3 The Label Approval Record (WI052-F1) will be given to a member of the Quality department for review.
- 5.4 A representative of the Maxx Orthopedics quality team will review label copy and compare to corresponding Production Procedures referenced in section 2.0 to ensure all necessary and correct information is on the labels.

WORK INSTRUCTION

TITLE:	Product Label Approval Process
---------------	---------------------------------------

- 5.5 The Maxx Orthopedics quality representative will also scan the QR codes to ensure the accuracy of information is contained in the QR codes as well.
- 5.6 If approved, the Maxx Orthopedics quality representative will sign in the “Quality Approval” section location at the bottom of the form.
- 5.7 The Quality Representative will then give the signed Label Approval Record back to the Production team for further processing.

Revision History

Revision	Details	ECO #	Author
01	Initial Release	25-068	D Angelichio